


## TO CHANGE ATTACHMENTS

2 custom massage head attachments available:

To change to a different set, simply twist off (counterclockwise) the set that is attached to the unit and twist on (clockwise) the desired set as shown.



## MAINTENANCE

### TO CLEAN:

Be sure to unplug the unit and allow it to cool before cleaning. Use a soft, slightly damp cloth to wipe. NEVER allow water or any liquids to come into contact with the handle.

### TO STORE:

Unplug the appliance from the outlet and allow it to cool before storing in its box or a clean, dry place.

DO NOT hang the unit by the power cord.


## SPECIFICATIONS

Name: Handhold Massager

Input: ☐ 220-240V ~ 50/60Hz ☐ 100-120V ~ 50/60Hz

Power: 26W

Note: specification are subject to change without notice.



503D-1 V.A



## HANDHOLD MASSAGER PROFESSIONAL PERCUSSION



Instruction Manual

## IMPORTANT SAFETY

**WHEN USING ELECTRICAL PRODUCTS, ESPECIALLY WHEN CHILDREN ARE PRESENT, BASIC SAFETY PRECAUTIONS SHOULD ALWAYS BE FOLLOWED, INCLUDING THE FOLLOWING:  
READ ALL INSTRUCTIONS BEFORE USING.**

**DANGER TO REDUCE THE RISK OF ELECTROCUTION:**

ALWAYS unplug the appliance from the electrical outlet immediately after using and before cleaning.

- DO NOT reach for an appliance that has fallen into water. Unplug it immediately.
- DO NOT use while bathing or showering.
- DO NOT place or store appliance where it can fall or be pulled into a tub.
- DO NOT place in or drop into water or other liquid.

**WARNING - TO REDUCE THE RISK OF BURNS, ELECTROCUTION, FIRE OR INJURY TO PERSONS:**

- An appliance should NEVER be left unattended when plugged in. Unplug from outlet when not in use and before putting on or taking off parts or attachments.
- Close supervision is necessary when this appliance is used by, on or near children, invalids or disabled persons.
- Do Not use outdoors.

- DO NOT use over sensitive skin areas or in the presence of poor circulation. Use this appliance only for its intended use as described in this manual.
- DO NOT use any attachments not provided with the unit.
- NEVER operate it if it has a damaged cord/plug, if it is not working properly, if it has been dropped into water or damaged.
- Keep cord away from heated surfaces.
- NEVER use while sleeping or fall asleep while using the massager.
- NEVER drop or insert any object into any opening.
- DO NOT operate where aerosol (spray) products are being used or where oxygen is being administered.
- DO NOT use on or near eyes or other highly sensitive areas.
- DO NOT carry this appliance by supply cord or use cord as handle.
- DO NOT use the percussion action massager directly on one concentrated area of the body. Percussion action massagers are designed to be used in a continuous motion, constantly changing the massage area.
- Turn all controls to the OFF position and remove plug from outlet.
- Choking hazard : small parts. Not recommended for children under 12 years of age.



**CAUTION:** —  
**pregnancy, menstruating or illness, consult your doctor before use.**



**CAUTION:**

**PLEASE READ ALL INSTRUCTIONS CAREFULLY BEFORE OPERATING.**

This product is not intended for medical use. Consult your doctor before using this product as below

- Pregnant
- with A pacemaker
- You have any concerns regarding your health
- Never leave the appliance unattended, especially if children are present.
- NEVER cover the appliance when it is in operation.
- DO NOT use this product for more than 15 minutes at a time.
- Extensive use could lead to the product's excessive heating and shorter life. Should this occur, discontinue use and allow the unit to cool before operating.
- NEVER use this product directly on swollen or inflamed areas or skin eruptions.
- DO NOT use it before bed. It has a stimulating effect and can delay sleep.
- NEVER use this product while in bed.
- This product should NEVER be used by any individual suffering from any physical ailment that would limit the user's capacity to operate the controls.
- This unit should not be used by children or invalids without adult supervision.


## How to use

**ALWAYS make sure the unit is set to the OFF position before inserting the plug into outlet.**

1. Massage: This product has two intensity settings (see Figure). To begin, push the massager intensity switch to the setting you want to use:

- a. Slow- For a gentle, soothing massage.
- b. Fast- Use for a more intense massage experience.

- 4. Apply massager head lightly to the area you wish to massage.
- 5. Move the unit slowly, holding on the targeted area. Massage the area for a short period of time and move to another area. NEVER massage the same spot longer than 3 minutes. CAUTION - as stated in the IMPORTANT SAFETY section of this manual, NEVER use this appliance on or near eyes or other highly sensitive areas.
- 6. Operation of longer than 15 minutes is not recommended and may cause overheating.
- 7. When massage is completed, turn unit off, unplug the power supply and store properly in a dry place.
- 8. Household use only.





Eases muscle aches and strain  
from your neck to your legs.



**OFF With Two Speed Settings** Choose from gentle or invigorating  
massage at the slide of a switch

**Infrared Ray Function** To soothe tired muscles

**Ergonomic Rubberized Handle**  
For easy, grip and handling



**Massage head**  
Tapping motion for a  
comforting massage



**2 Custom Attachments**  
Each provides a unique  
massage sensation

All right reserved, any alterment please prefer to the real product.