

IMPORTANT SAFEGUARDS

**WHEN USING ELECTRICAL PRODUCTS, ESPECIALLY WHEN CHILDREN ARE PRESENT, BASIC SAFETY PRECAUTIONS SHOULD ALWAYS BE FOLLOWED, INCLUDING THE FOLLOWING:
READ ALL INSTRUCTIONS BEFORE USING
DANGER- TO REDUCE THE RISK OF ELECTROCUTION:**

- Always unplug the appliance from the electrical outlet immediately after using and before cleaning.
- DO NOT reach for an appliance that has fallen into water. Unplug it immediately.
- DO NOT use while bathing or in shower.
- DO NOT place in or drop into water or other liquid.
- NEVER use pins or other metallic fasteners with this appliance.
- Keep dry - DO NOT operate in a wet or moist condition.
- DO NOT place or store appliance where it can fall or be pulled into a tub or sink.

WARNING - TO REDUCE THE RISK OF BURNS, ELECTROCUTION, FIRE OR INJURY TO PERSONS:

- An appliance should never be left unattended when plugged in. Unplug from outlet when not in use and before putting on or taking off parts or attachments.
- Close supervision is necessary when this appliance is used by, on or near children, invalids or disabled persons.
- Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by our company; specifically any attachments not provided with the unit.
- NEVER operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water.
- Keep cord away from heated surfaces.
- NEVER use while sleeping or fall asleep while using the massager.
- NEVER drop or insert any object into any opening.
- DO NOT operate where aerosol (spray) products are being used or where oxygen is being administered.
- DO NOT operate under a blanket or pillow. Excessive heating can occur and cause fire, electrocution or injury to persons.
- DO NOT carry this appliance by the power cord or use cord as handle.
- DO NOT use on an infant or invalid or on a sleeping.

- DO NOT crush. Avoid sharp folds.
- DO NOT use outdoors.
- DO NOT use on insensitive skin or on a person with poor blood circulation.

SAVE THESE INSTRUCTIONS

Notes

- The operation will be stopped automatically by the auto-timer function after 10 minutes, do not forget to switch off the Main Switch after leaving.
- After using the unit for 40 minutes continuously, allow a break for massager to prevent overheating.



Please read all instructions carefully before operating.

- It is intended only to provide a luxurious massage.
- Consult your doctor prior to using this product, if
 - pregnant
 - with a pacemaker
 - You have any concerns regarding your health
- Diabetics.
- Never use this product directly on swollen or inflamed areas.
- DO NOT use this product before or in bed. It has a stimulating effect and can delay sleep.
- This product should NEVER be used by any individual suffering from any physical ailment that would limit the user's capacity to operate the controls or who has sensory deficiencies in the lower half of their body.

MAINTENANCE

To Store

Place massager in its box or in a safe, dry, cool place. Avoid contact with sharp edges or pointed objects which might cut or puncture the fabric surface. To avoid breakage, do not wrap the power cord around the unit. Do not hang the unit by the hand control cord.

To Clean

Unplug the unit and allow it to cool before cleaning. Spot clean only with a soft, slightly damp sponge. Never allow water or any other liquids to come into contact with the unit.

- Do not immerse in any liquid to clean.
- Never use abrasive cleaners, brushes, gasoline, kerosene, glass / furniture polish or paint thinner to clean.

Benefits:

- Effective at relaxing shoulder / neck muscles and easing aches and pains.
- Suitable for mid/ lower back, abdominal, calf and thigh massage.
- Improve blood circulation and lymphatic flow.
- Light and portable massage wherever you are, at home or in the office.
- Reduce stress.




How To Use


Insert the power plug into electric outlet.

- Drape it across your shoulders and you will be targeting your stiff shoulders.
- Drop it back slightly and you can feel the tapping massage reaching your mid back.
- Lie down on your tummy and feel the drumming motion delivering comfort to your lower back.
- Wrap it around your thighs, calves or knees and experience the benefits of improved blood circulation to under-exercised regions of your body.
- Heat therapy: It is effective at relieving aches and pain.
- Intensity Level: 8 different intensity levels available.

(+) UP (-) DOWN



Control Panel



A TRULY HANDS-FREE TAPPING MASSAGER

Relieving aches and pains with heat therapy.

Instruction Manual

Key Features:

- Easy use with control panel
- 8 different intensity levels
- Soothing heat therapy
- 41 Unique rhythmic massage cycles

SPECIFICATIONS

Name: HANDS-FREE TAPPING MASSAGER

Model: DF703

Input: □ AC220-240V ~ 50/60Hz

□ AC100-120V ~ 50/60Hz

Power: 60W

